




# Duramark® Wristbands



SATO Life Science ID Solutions



Infant Wristband



Adult Wristband

## SATO's DuraMark® Wristbands Enhanced Safety Through Positive Identification

SATO's DuraMark® Wristbands are virtually indestructible and display various details of patient information, encoded securely in a barcode. They remain legible even after extended exposure to water, abrasion, cleaners and everyday wear.


Each wristband provides enhanced safety through positive identification of individual patients, linked to a patient's medical or account number using both human-readable text, barcodes and if desired, a patient's picture. All without the need of thermal transfer ribbons which raise costs, are cumbersome, and under HIPAA guidelines require expensive destruction procedures.

### SATO's DuraMark® Wristbands

- Smudge, Smear, and Water Proof
- Comfortable & Dependable
- Latex Free Material

### SATO's Failsafe Fasteners™

- 7 Colors For Identity, Status Codes
- Non-Transferable Security Closure
- Latex Free Plastic





### Wristband Product Information

These all-purpose bands provide reliable, sharp-edge print quality to maximize the barcode scan accuracy rate for patient identification.

DESCRIPTION	ITEM #	SIZE	COLOR
<b>Adult Wristband</b>	60S101002	Starting Size: 1.535" W x 12" L Finished Size: 1.417" W x 11.575" L	Bright White
<b>Infant Wristband</b>	60S121001	Starting Size: 1.545" W x 7.75" L Finished Size: .75" W x 7.35" L	Bright White

PROPERTIES	
<b>Construction</b>	Multi-layered non-cavitated films
<b>Sensing</b>	Eyemark
<b>Imaging Type</b>	Direct Thermal
<b>Scannability</b>	Visible, NIR and IR light spectrums
<b>FDA Approvals</b>	Complies with 21CFR175.105
DURABILITY	
<b>Smudge &amp; Scratch Test</b>	Exceeds 1000 rubs using crockmeter
<b>UV Light Resistance</b>	Minimal discoloration with extended UV exposure
<b>Water Resistance</b>	Withstands total emersion and wet heat up to 150°F (66°C)
<b>Heat Resistance</b>	1 hour at 185°F (85°C)
<b>Alcohol Resistance</b>	No reaction to alcohol
<b>Storage Conditions</b>	41 - 77°F (5 - 25°C) RH 45-65%
<b>Storage Life</b>	1 Year after manufacturing
<b>Usage Conditions</b>	41 - 104°F (5 - 40°C)
PACKAGING	
<b>Roll Size</b>	O.D. 4.4" / I.D. 1.5 inches with fiber core
<b>Roll Count</b>	100 wristbands per roll
<b>Roll Wind</b>	Face Out
<b>Rolls per Box</b>	5 rolls
<b>Clips Sold Separately</b>	520 clips per box, per color Available in White, Green, Yellow, Pink, Purple, Blue and Red

The above data represents product averages, allowing for industry-accepted variances. These products should be tested in the end-use conditions to ensure that they meet the requirements of the specific application.



WS4 DT  
Printer



**CG and WS4 Series Printers  
are Cerner® Validated**



CG2 DT  
Printer



**SATO AMERICA  
CORPORATE**  
10350-A Nations Ford Rd  
Charlotte, NC 28273  
Phone: (704) 644-1650  
Fax: (704) 644-1662  
sales-sallc@sato-global.com

**LABEL MANUFACTURING, SERVICE & SALES**  
UNITED STATES MEXICO CENTRAL AMERICA  
CANADA CARIBBEAN



[www.satoamerica.com/healthcare](http://www.satoamerica.com/healthcare)